



## Oyster Tabmaster – Bolted Body

**The flexible leaf is co-moulded in FDA compatible silicone and gives the flexibility necessary to ensure no damage to tablets or capsule being processed.**

All Oyster™ branded sanitary valves are designed according to cGMP standards, for trouble free operation and easy maintenance. The valve design ensures optimal hygienic and functional characteristics and eliminates recesses with hard-to-clean corners. Oyster™ branded valves are manufactured to exacting tolerances and finished to the highest pharmaceutical standards.

APPLICATION TYPE



Shut-off



## FEATURES

Atex Certification N.A.

Diameter DN 100 / 150 / 200 / 250 – (4”-6”-8”-10”)

WIP suitable for washing in place

Operation Hand lever  
Pneumatic actuator

## SPECIFICATIONS

M.O.C. CONTACT PARTS AISI 316L, C22, (other alloys available upon request)

M.O.C. OTHER PARTS AISI 304

M.O.C. LEAF Co-moulded Silicone (FDA 177.2600 compliant)

M.O.C. SEALS (all FDA  
177.2600 compliant) SILICONE, EPDM

AVAILABLE CONNECTIONS TC BS 4825-3, TC ASME BPE, TC DIN32676, TC ISO2852, TC  
OSAME A, Weld End, Connecting Flange, Wafer